

THE NEW
CapnoPen[®]



Article no. CP-001k

The only **PIPAC[®]** device certified for aerosolization of chemotherapy, DNA, RNA, nanomolecules, viruses and immunotherapy into body cavities

Made by [CapnoPharm](#) (legal manufacturer)

CE certification (MDR)


CapnoPharm GmbH

Vor dem Kreuzberg 17 . D-72070 Tübingen . Deutschland (Germany)


+49 (0) 7071 / 855 07 45 | +49 (0) 7071 / 855 07 46 | info@capnopharm.com | www.capnopharm.com

Optimal driving pressure, spray angle and droplet size

CapnoPen[®] CP-001k



Previous CapnoPen[®]



Superior
spray angle

∠ 57° to 66°

Optimal
droplet size

💧 (29.1 ± 1.5µm)

Preserves the peritoneum

- no serosal injuries
- includes a protective filter
- no metal oxide residue in the target tissue

High-quality manufacturing

- German engineering and manufacturing
- only approved material (100% stainless steel)

Safe design

- smooth, high-end surface
- conical safety tip
- prevents tubing disconnection
- high-end sterile packaging

100% quality control

- quality control using artificial intelligence
- each CapnoPen[®] tested for spray pattern and tightness

Competitive pricing

Technical support in 40+ countries



Scan for more information
www.capnopharm.com

The product is not registered and approved for sale in all countries. Indication for use may also vary depending on the country. All images in this brochure are for reference only. The product may only be used by trained physicians. For further information and technical details please contact us. Subject to technical changes without further notice.

CapnoPharm GmbH

Vor dem Kreuzberg 17 . D-72070 Tübingen . Deutschland (Germany)

☎ +49 (0) 7071 / 855 07 45 | ☎ +49 (0) 7071 / 855 07 46 | ✉ info@capnopharm.com | 🌐 www.capnopharm.com