



38800Q

NON-MAGNETIC
MRI PATIENT MONITOR



 **IRADIMED**

THE IRADIMED 3880Q

WORLD'S ONLY NON-MAGNETIC MRI MONITOR

This one of a kind, lightweight, portable monitor has been meticulously engineered to meet the needs of today's complex MRI workflow. **From bedside through transport, the 3880Q non-magnetic monitor is mountable patient-side, on a wall, roll-stand, bed rail, or an anesthesia cart.** The 3880Q can be quickly detached for immediate mobility and provides MRI safety with full functionality in a compact, non-magnetic package.



NON-MAGNETIC DESIGN

- 30,000 Gauss Field Line safety
- Zone IV piece of mind
- Eliminates projectile risks



IMPROVED WORKFLOW

- 10 lbs (4.5kg), small, lightweight design
- Familiar bedside monitor look and feel
- Flexible mounting: cart, IV pole & bed



CONTINUOUS MONITORING

- Interdepartmental transfer ease
- Uninterrupted vital sign monitoring
- Patient safety is enhanced

ROADMAP TO A SEAMLESS WORKFLOW



1. 3880Q PATIENT SETUP

Connecting the 3880Q MRI Patient Monitor to the patient within the 'safety-net' of their care unit **improves patient stability** prior to their MRI appointment.




2. PATIENT TRANSPORT TO MRI

The lightweight design of the 33880Q MRI Patient Monitor allows minimal members to **easily transport the patient to the MRI** without the need to transfer monitors again upon arrival.



4. RECOVERY OR CRITICAL CARE

Whether the patient's post MRI transfer is to their originating department or a recovery area, the 3880Q MRI Patient Monitor is with the patient every step of the way, **standardizing continuity of care.**



3. MRI EXAM

With its 30,000 gauss rating and small footprint, clinicians have the freedom to ergonomically position the 3880Q MRI Patient Monitor **where it best enhances patient care** and throughput.

MODERN PROBLEMS REQUIRE MODERN SOLUTIONS



WIRELESS ECG WITH DYNAMIC GRADIENT FILTERING

- Industry exclusive "Convertible ECG" supports 3, 4, and 5 leads setups giving clinicians flexibility on patient application.
- Single lead wire replacement instead of whole cable reduces replacement cost.



WIRELESS MASIMO® PULSE OXIMETRY

- More than 100 studies have shown that Masimo® SET provides best in class performance with patients with low perfusion.
- Integrated motion tolerance provides stability during transport results in more reliable readings.



WIRELESS INVASIVE BLOOD PRESSURE MONITORING

- Intuitive design and setup keeps lines on the patient's bed and off the floor.
- Reduces the risks of tripping hazards on long IBP lines.




NIBP WITH ecQ-TEK

- Quick inflation cycles.
- Energy efficient ecQ-TEK pump technology extends battery run time.
- Quiet and smooth operation enhances reliability and comfort.



INTEGRATED BATTERY CHARGING

- No external Zone III battery chargers required to manage.
- Assures wireless PODS always have a charge and are ready for immediate use.



UNINTERRUPTED PATIENT TEMPERATURE

- Fiber optic thermometer provides real time and reliable readings.
- Precision results within 0.3° C instills confidence during clinical decisions.



LOW FLOW END TIDAL CO₂


- Magnetic field doesn't alter accuracy or reliability giving clinicians peace of mind.
- Cross contamination is controlled by the enclosed design, eliminating the need for a water-trap.



WIRELESS ANESTHESIA GAS MONITORING


- Operating Room workflow supported by mounting on the anesthesia machine.
- Multiple 3880Q patient monitors supported by a single multi-gas module.

PERFECTLY PORTABLE: WIRELESS MONITORING BY REMOTE CONTROL



EXTENDED RANGE WIRELESS REMOTE

- Remotely monitor patients inside the MRI Zones III & IV with minimal interruptions.
- Tablet wirelessly controls the 3880Q with full audio and visual alarm functionality.



E.M.R CONNECTIVITY THAT SAVES TIME

- HL7 and RS232 simultaneous output options allows integration to both legacy and future EMR systems.

TOUCH SCREEN TECHNOLOGY

- Every step of the patient's MRI care cycle has been simplified with an intuitive touch screen user interface.




SECURE 4K REMOTE VIEWING

- The Remote Viewing Station allows for clinicians to oversee up to four patients on a single screen.



Reverse side of 3885B



Remote monitoring 4 patients on 1 screen

SERVICE BEYOND REPAIR!

DEPOT SERVICE GETS YOU BACK TO PATIENT CARE FAST!

IRadimed depot service combined with our priority equipment loaner program is the most effective way to ensure MRI cases are resumed promptly when unexpected service is needed.



No one knows MRI patient care devices like IRadimed. IRadimed technical support and medical device service engineers are expert problem solvers with deep product knowledge. Our dedicated team of professionals are eager to manage the end-to-end care of your 3880Q MRI Patient Monitor System and accessories.

Scan for Training Videos



www.IRadimed.com



UPGRADE YOUR PATIENT CARE BY GOING NON-MAGNETIC!

Start moving forward by replacing your magnetic, bulky, outdated equipment with a portable, lightweight, MRI monitor. From patient preparation to recovery, the 3880Q MRI patient-side monitor provides your facility with more reliability, and functionality, in a compact, non-magnetic package.

*Not all features are available in all markets.

 **IRADIMED**

1025 Willa Springs Drive
Winter Springs, FL 32708
(407) 677-8022

