













IRadimed 3880 MRI Patient Monitoring System

Accessory Guide

ECG Accessories



Part Number	Description	Image
1812-25	Hexrode Electrode (Box of 25) <i>MRI electrode for 3, 4 and 5 lead applications.</i>	
9009	ECG Skin Prep Gel (Box of 3)	
1811	9" Five Lead ePOD Leadwire (AHA/AAMI) <i>ECG leadwire for 3 and 5 lead ECG applications with AAMI color coding.</i>	
1811E	9" Five Lead ePOD Leadwire (IEC) <i>ECG leadwire for 3 and 5 lead ECG applications with IEC color coding.</i>	
3881	8 Channel Wireless EPOD, MRI ECG	

SpO2 Accessories








Part Number	Description	Image
1822A-25	Adult SpO2 Secure Wrap Grip Sensor (Box of 25) <i>Disposable SpO2 applicator designed for adult digits.</i>	
1822P-25	Pediatric SpO2 Secure Wrap Grip Sensor (Box of 25) <i>Disposable SpO2 applicator designed for pediatric digits.</i>	
1822I-25	Infant SpO2 Secure Wrap Grip Sensor (Box of 25) <i>Disposable SpO2 applicator designed for infant digits.</i>	
1822N-25	Neonatal SpO2 Secure Wrap Grip Sensor (Box of 25) <i>Disposable SpO2 applicator designed for neonatal hands.</i>	
1171	Reusable Grip Sensor (Size 1, 2 and 3)	
1821	Fiber Optic oPOD SpO2 Cable	
3882	8 Channel Wireless oPOD, SpO2 Modul incl. Grip Sensor Kit (1171) & SpO2 Cable (1821)	




NIBP Accessories

Part Number	Description	Image
1832LA	Large Adult Reusable NiBP Cuff (34 - 44 cm)	
1832A	Adult Reusable NiBP Cuff (27 - 35 cm)	
1832SA	Small Adult Reusable NiBP Cuff (20.5 - 28.5 cm)	
1832P	Pediatric Reusable NiBP Cuff (14 - 21.5 cm)	
1832I	Infant Reusable NiBP Cuff (9 - 15 cm)	
1833N4-10	Neonatal Size 4 NiBP Cuff (Box of 10) (7 - 13 cm) <i>Deliverytime: 2 weeks</i>	
1833N3-10	Neonatal Size 3 NiBP Cuff (Box of 10) (6 - 11 cm) <i>Deliverytime: 2 weeks</i>	
1833N2-10	Neonatal Size 2 NiBP Cuff (Box of 10) (4 - 8 cm) <i>Deliverytime: 2 weeks</i>	

1833N1-10	Neonatal Size 1 NiBP Cuff (Box of 10) (3 - 6 cm) <i>Deliverytime: 2 weeks</i>	
1831	Blood Pressure Hose (NIBP)	

Respiration Accessories




Part Number	Description	Image
1842A-25	Adult Nasal Cannula (Pack of 25)	
1842P-25	Pediatric Nasal Cannula (Pack of 25)	
1842I-25	Infant Nasal Cannula (Pack of 25)	
1843A-25	Adult Oral / Nasal Cannula (Pack of 25)	
1843P-25	Pediatric Oral / Nasal Cannula (Pack of 25)	
1844A-25	Adult Divided Nasal Cannula (Pack of 25)	
1844P-25	Pediatric Divided Nasal Cannula (Pack of 25)	


1841-25	Filtered Sample Line Extension, 2.4m (Pack of 25) <i>Coextruded sample line with hydrophobic filter. Deliverytime: 2 weeks</i>	
1849-25	Nomoline Agent Sample Line, 2.0m (Box of 25) <i>For use with IRadimed's 3886 Wireless Multigas Module only. Has a male luer lock. 1841 sample lines can be used to extend this when needed. Deliverytime: 2 weeks</i>	
1850-25	Nomoline Agent Reusable Filter (Box of 25) <i>For use with IRadimed's 3886 Wireless Multigas Module only. Must be used with 1841 IRadimed sample lines. Deliverytime: 2 weeks</i>	

Temperature Accessories





Part Number	Description	Image
1851	Optical Temperature Sensor	

Mounting Accessories

Part Number	Description	Image
1891	3880 Flexible Mounting Adapter	
1119	IRadimed Roll Stand	
1883	3886 Mounting Arm	

1884	3880 Mounting Arm	
-------------	--------------------------	---

Miscellaneous Accessories

Part Number	Description	Image
1882-3	3885B Recorder Paper (Pack of 3)	
1881	Universal Gating Interface	
1133	3880 Battery <i>Non-Magnetic battery used with the 3880</i>	
1188	3885T Battery <i>Non-Magnetic battery used with the 3885T</i>	
1887	Dual Radio Module <i>Allows an additional 3880 system to communicate with an IRadimed external add-on accessory such as Wireless Multi-Gas Module and/or IBP Module.</i>	