

CONSUMABLES

Single and double head injectors for computed tomography (21 bar/305 psi)



Article number

ELS 200 ml (Easy Loading Syringe) with filling tube 31 6025-000

Residual volume: 1.5 ml

ELS 200 ml (S) (Easy Loading Syringe) 31 6026-000

Residual volume: 1.5 ml

Packing unit: 50 pieces/carton – 1,500 pieces/pallet



ES 223 tube system with spike 31 7080-100

	ID	Length	Filling volume
Pressure side	2.0 mm	13 cm	0.4 ml
Suction side	2.7 mm	100 cm	6.3 ml

Packing unit: 50 pieces/carton – 2,800 pieces/pallet

just for
Accutron CT



EK 223 tube system with drip chamber 31 7081-100

	ID	Length	Filling volume
Pressure side	2.0 mm	13 cm	0.4 ml
Suction side	2.7 mm	100 cm	5.8 ml
Drip chamber			10 ml

Packing unit: 50 pieces/carton – 2,800 pieces/pallet

just for
Accutron CT



Tube system with suction line and drip chamber (for Scanbag® CM-supply bag) 31 7083-100

	ID	Length	Filling volume
CM-Side	2.0 mm	10 cm	0.3 ml
NaCl-Side	2.0 mm	40 cm	1.3 ml
Patient line	1.5 mm	1.500 mm	2.7 ml
Suction side with drip chamber	2.7 mm	1.000 mm	15.8 ml

Suitable for prefilled syringes

Packing unit: 50 pieces/carton – 2,800 pieces/pallet

CONSUMABLES

Single and double head injectors for computed tomography (21 bar/305 psi)



MRS 222 MR-tube system 31 7099-100 with 1 inline drip chamber (for Scanbag® CM-supply bag)

Pressure side, ID 2.0 mm	Length	Filling volume
CM-Side	12 cm	0.6 ml
NaCl-Side	36 cm	1.3 ml
Suction side, ID 2.7 mm	Length	Filling volume
CM-Side, green	100 cm	5.8 ml
NaCl-Side, white	100 cm	5.8 ml
2 x drip chamber		each one 10 ml

Packing unit: 50 pieces/carton – 2,800 pieces/pallet



MRS 222 MR-tube system 31 7100-100

Pressure side, ID 2.0 mm	Length	Filling volume
CM-Side	12 cm	0.6 ml
NaCl-Side	36 cm	1.3 ml
Suction side, ID 2.7 mm	Length	Filling volume
CM-Side, green	100 cm	5.8 ml
NaCl-Side, white	100 cm	5.8 ml
2 x drip chamber		each one 10 ml

Packing unit: 50 pieces/carton – 2,800 pieces/pallet



Double head-tube system 31 7101-100

Length from Y-Connector: 180 cm
Inner diameter: 2.0 mm
Complete filling volume: 6.9 ml

Packing unit: 50 pieces/carton – 4,000 pieces/pallet



Double head tube system with 31 7109-100 1 inline drip chamber (for Scanbag® CM-supply bag)

Pressure side, ID 2.0 mm	Length	Filling volume
CM-Side	12 cm	0.4 ml
NaCl-Side	36 cm	1.1 ml
Patient line (ID: 1.5 mm)	150 cm	2.6 ml
Suction side, ID 2.7 mm	Length	Filling volume
CM-Side, green	100 cm	5.8 ml
NaCl-Side, white	100 cm	5.8 ml
2 x drip chamber		each one 10 ml

Packing unit: 50 pieces/carton – 2,800 pieces/pallet

CONSUMABLES

Single and double head injectors for computed tomography (21 bar/305 psi)



MRS 225 MR-tube system

31 7125-100

Inner diameter: 2.0 mm

	Length	Filling volume
CM-Side	6 cm	0.2 ml
NaCl-Side	32 cm	1.0 ml

Packing unit: 50 pieces/carton – 4,000 pieces/pallet



82 CT Kit ELS

31 7616-000

1 piece ELS 200 ml (Easy Loading Syringe) with filling tube 31 6025-000
1 piece Coiled line 150 cm 31 8213-000

Packing unit: 50 pieces/carton – 1,500 pieces/pallet

just for
Accutron CT



CT2/MRT Set ELS

31 7625-000

2 pieces ELS 200 ml (Easy Loading Syringe) with filling tube 31 6025-000
1 piece MRS 222 MR-tube system 31 7100-000
1 piece ES 224/150 patient line 31 8151-000

Packing unit: 20 pieces/carton – 600 pieces/pallet



CT2/MRT Set ELS (S)

31 7627-100

2 pieces ELS 200 ml (S) (Easy Loading Syringe) 31 6026-000
1 piece MRS 222 MR-tube system 31 7100-100
1 piece ES 224/150 patient line 31 8151-100

Packing unit: 20 pieces/carton – 600 pieces/pallet

CONSUMABLES

Single and double head injectors for computed tomography (21 bar/305 psi)

Article number



ES 224/25 patient line with valve

31 8020-100

Length: 25 cm
Innen diameter: 1.5 mm
Filling volume: 0.4 ml

Packing unit: 100 pieces/carton – 8,000 pieces/pallet



One-way-stopcock with
extension line 30cm, luer lock

31 8025-100

Length: 32 cm
Inner diameter: 2.0 mm
Filling volume: 1.2 ml

Packing unit: 200 pieces/carton – 5,000 pieces/pallet



Patient line 150 cm
with terminal check valves

31 8026-100

Length: 150 cm
Inner diameter: 2.0 mm
Filling volume: 4.8 ml

Packing unit: 50 pieces/carton – 5,850 pieces/pallet



Patient line 150 cm with 2 check valves

31 8036-100

Length: 150 cm + 10 cm
Inner diameter: 2.0 mm
Filling volume: 5.2 ml

Packing unit: 50 pieces/carton – 5,850 pieces/pallet

CONSUMABLES

Single and double head injectors for computed tomography (21 bar/305 psi)

Article number



DW 229 Connector with 3-way stop cock

31 8042-100

Filling volume: 0.5 ml

Packing unit: 100 pieces/carton – 8,000 pieces/pallet



ES 224 Patient line Inner diameter: 1.5 mm

ES224/80, length: 80 cm, filling volume: 1.4 ml

31 8081-100

ES224/100, length: 100 cm, filling volume: 1.8 ml

31 8101-100

ES224/150, length: 150 cm, filling volume: 2.7 ml

31 8151-100

ES224/200, length: 200 cm, filling volume: 3.6 ml

31 8199-100

ES224/250, length: 250 cm, filling volume: 4.4 ml

31 8250-100

ES224/300, length: 300 cm, filling volume: 5.3 ml

31 8300-100

Packing unit: 200 pieces/carton – 11,200 pieces/pallet



Patient line with valve

Inner diameter: 1.5 mm

Length: 150 cm, filling volume: 2.7 ml

31 8152-100

Length: 200 cm, filling volume: 3.5 ml

31 8198-100

Length: 250 cm, filling volume: 4.4 ml

31 8251-100

Length: 300 cm, filling volume: 5.3 ml

31 8301-100

Packing unit: 200 pieces/carton – 11,200 pieces/pallet



Coiled line 150 cm

31 8213-000

Length: 150 cm

Inner diameter: 1.6 mm

Filling volume: 3.1 ml

Packing unit: 50 pieces/carton – 5,850 pieces/pallet

CONSUMABLES

Single and double head injectors for computed tomography (21 bar/305 psi)

Article number



Coiled line with check valve

Inner diameter: 1.6 mm


Length: 150 cm, filling volume: 3.1 ml

Length: 180 cm, filling volume: 3.7 ml

31 8214-000

31 8181-000

Packing unit: 50 pieces/carton - 5,850 pieces/pallet



Coiled line 150 cm with 2 check valves


31 8217-000

Inner diameter: 1.5 mm

Length: 150 cm

Filling volume: 2.6 ml

Packing unit: 50 pieces/carton - 5,850 pieces/pallet



Coiled line 235 cm with Y-connector and valve

31 8235-000

Length: 235 cm + 1 x 10 cm + 1 x 12 cm

Inner diameter: 1.6 mm

Total Filling volume: 5.6 ml

Packing unit: 50 pieces/carton - 5,850 pieces/pallet



Coiled line with check valve

Inner diameter: 1.5 mm

Length: 250 cm, filling volume: 4.4 ml

Length: 300 cm, filling volume: 5.3 ml

31 8253-000

31 8303-000

Packing unit: 50 pieces/carton - 5,850 pieces/pallet

CONSUMABLES

Single and double head injectors for computed tomography (21 bar/305 psi)



SD 225 suction tube with 3-way stopcock and spike

Article number

31 9080-100

Filling volume: 0.3 ml (3-way stopcock) + 6.3 ml (suction tube)

Packing unit: 50 pieces/carton – 4,000 pieces/pallet

just for
Accutron CT



DW 226 3-way stopcock

31 9085-100

Filling volume: 0.1 ml

Packing unit: 50 pieces/carton
200 pieces/carton – 11,200 pieces/pallet



SP 227 vented spike, luer lock

31 9095-100

Filling volume: 0.4 ml

Packing unit: 50 pieces/carton
200 pieces/carton – 11,200 pieces/pallet



Minispike, green

31 9099-100

Filling volume: 0.2 ml

Packing unit: 50 pieces/carton
200 pieces/carton – 11,200 pieces/pallet