

CONSUMABLES

Double head injectors for magnetic resonance imaging (21 bar/305 psi)



MRS 222 MR-tube system

Article number

31 7100-100

Pressure side, ID 2.0 mm	Length	Filling volume
CM-Side	12 cm	0.6 ml
NaCl-Side	36 cm	1.3 ml
Suction side, ID 2.7 mm	Length	Filling volume
CM-Side, green	100 cm	5.8 ml
NaCl-Side, white	100 cm	5.8 ml
2 x drip chamber		each one 10 ml

Packing unit: 50 pieces/carton – 2,800 pieces/pallet



Double head-tube system

31 7101-100

Length from Y-piece: 180 cm
 Inner diameter: 2.0 mm
 Total filling volume: 6.9 ml

Packing unit: 50 pieces/carton – 4,000 pieces/pallet



MR-ELS tube system with 2 valves, green

31 7104-000

Inner diameter: 1.7 mm
 Coiled line: 180 cm
 Short side: 10 cm
 Long side: 34 cm
 Filling volume: 5.1 ml

Packing unit: 50 pieces/carton



Double head-tube system with 1 inline drip chamber

31 7109-100

(for Scanbag® CM-supply bag)

Pressure side, ID 2.0 mm	Length	Filling volume
CM-Side	12 cm	0.4 ml
NaCl-Side	36 cm	1.1 ml
Patient line (ID: 1.5 mm)	150 cm	2.6 ml
Suction side, ID 2.7 mm	Length	Filling volume
CM-Side, green	100 cm	5.8 ml
NaCl-Side, white	100 cm	5.8 ml
2 x drip chamber		each one 10 ml

Packing unit: 50 pieces/carton – 2,800 pieces/pallet

Double head injectors for magnetic resonance imaging (21 bar/305 psi)



MRS 225 MR-tube system


Article number

31 7125-100

Inner diameter: 2.0 mm

	Length	Filling volume
CM-Side	6 cm	0.2 ml
NaCl-Side	32 cm	1.0 ml

Packing unit: 50 pieces/carton – 4,000 pieces/pallet



MR-Set ELS 65 ml

31 7602-100

2 pieces	ELS 65 ml (Easy Loading Syringe)	31 6065-000
1 piece	MR-ELS tube system with 2 valves, green	31 7104-100
2 pieces	standard spike	31 9095-100

Packing unit: 20 pieces/carton

Suitable for 65/65 ml version




MR-Set ELS 65/200 ml

31 7603-100

1 piece	ELS 65 ml (Easy Loading Syringe)	31 6065-000
1 piece	ELS 200 ml (S) (Easy Loading Syringe)	31 6026-000
1 piece	Double head tube system	31 7101-100
1 piece	SP 227 vented spike, luer lock	31 9095-100
1 piece	Minispike, green	31 9099-100

Packing unit: 24 pieces/carton – 384 pieces/pallet

Suitable for 65/200 ml version



MR-Kit ELS 65 ml

31 7604-100

1 piece	ELS 65 ml syringe	31 6065-000
1 piece	MR-ELS tube system with 2 valves, green	31 7104-100
1 piece	standard spike	31 9095-100

Packing unit: 50 pieces/carton

CONSUMABLES

Double head injectors for magnetic resonance imaging (21 bar/305 psi)



CT2/MRT Set ELS

Article number
31 7625-000

2 pieces	ELS 200 ml (Easy Loading Syringe) with filling tube	31 6025-000
1 piece	MRS 222 MR-tube system	31 7100-000
1 piece	ES 224/150 patient line	31 8151-000

Packing unit: 20 pieces/carton – 600 pieces/pallet

Suitable for 200/200 ml version

just for
Accutron MR



CT2/MRT Set ELS (S)

Article number
31 7627-100

2 pieces	ELS 200 ml (S) (Easy Loading Syringe)	31 6026-000
1 piece	MR 222 MR-tube system	31 7100-100
1 piece	ES 224/150 patient line	31 8151-100

Packing unit: 20 pieces/carton – 600 pieces/pallet

Suitable for 200/200 ml version

just for
Accutron MR



ES 224/25 patient line with valve

Article number
31 8020-100

Length: 25 cm
Inner diameter: 1.5 mm
Filling volume: 0.4 ml

Packing unit: 100 pieces/carton – 8,000 pieces/pallet



One-way-stopcock with extension line 30cm, luer lock

Article number
31 8025-100

Length: 32 cm
Inner diameter: 2.0 mm
Filling volume: 0.9 ml

Packing unit: 200 pieces/carton – 5,000 pieces/pallet

Double head injectors for magnetic resonance imaging (21 bar/305 psi)

Article number



Patient line 150 cm with terminal check valves

31 8026-100

Length: 150 cm
Inner diameter: 2.0 mm
Filling volume: 4.8 ml

Packing unit: 50 pieces/carton – 5,850 pieces/pallet



Patient line 150 cm with 2 check valves

31 8036-100

Length: 150 cm + 10 cm
Inner diameter: 1.5 mm
Filling volume: 2.8 ml

Packing unit: 50 pieces/carton – 5,850 pieces/pallet



DW 229 Connector with 3-way stop cock

31 8042-100

Filling volume: 0.2 ml

Packing unit: 100 pieces/carton – 8,000 pieces/pallet



ES 224 Patient line

Inner diameter: 1.5 mm

ES 224/80 Length: 80 cm, filling volume: 1.5 ml 31 8081-100

ES 224/100 Length: 100 cm, filling volume: 1.8 ml 31 8101-100

ES 224/150 Length: 150 cm, filling volume: 2.7 ml 31 8151-100

ES 224/200 Length: 200 cm, filling volume: 3.6 ml 31 8199-100

ES 224/250 Length: 250 cm, filling volume: 4.4 ml 31 8250-100

ES 224/300 Length: 300 cm, filling volume: 5.3 ml 31 8300-100

Packing unit: 200 pieces/carton – 11,200 pieces/pallet

CONSUMABLES

Double head injectors for magnetic resonance imaging (21 bar/305 psi)

Article number



Patient line with valve

Inner diameter: 1.5 mm

Length: 150 cm, filling volume: 2.7 ml

Length: 200 cm, filling volume: 3.5 ml

Length: 250 cm, filling volume: 4.4 ml

Length: 300 cm, filling volume: 5.3 ml

31 8152-100

31 8198-100

31 8251-100

31 8301-100

Packing unit: 200 pieces/carton – 11,200 pieces/pallet



Coiled line 150 cm

31 8213-000

Length: 150 cm

Inner diameter: 1.6 mm

Filling volume: 3.1 ml

Packing unit: 50 pieces/carton – 5,850 pieces/pallet



Coiled line with check valve

Inner diameter: 1.6 mm


Length: 150 cm, filling volume: 3.1 ml

Length: 180 cm, filling volume: 3.7 ml

31 8214-000

31 8181-000

Packing unit: 50 pieces/carton – 5,850 pieces/pallet



Coiled line 150 cm with 2 check valves

31 8217-000

Inner diameter: 1.5 mm

Length: 150 cm


Filling volume: 2.6 ml

Packing unit: 50 pieces/carton – 5,850 pieces/pallet

CONSUMABLES

Double head injectors for magnetic resonance imaging (21 bar/305 psi)

Article number



Coiled line 235 cm with Y-connector and valve 31 8235-000

Length: 235 cm + 1 x 10 cm + 1 x 12 cm

Inner diameter: 1.6 mm

Total Filling volume: 5.6 ml

Packing unit: 50 pieces/carton – 5,850 pieces/pallet

just for
Accutron MR3



Coiled line with check valve

Inner diameter: 1.5 mm

Length: 250 cm, filling volume: 4.4 ml

31 8253-000

Length: 300 cm, filling volume: 5.3 ml

31 8303-000

Packing unit: 50 pieces/carton – 5,850 pieces/pallet



DW 226 3-way stopcock

31 9085-100

Filling volume: 0.3 ml

Packing unit: 50 pieces/carton

200 pieces/carton – 11,200 pieces/pallet

CONSUMABLES

Double head injectors for magnetic resonance imaging (21 bar/305 psi)



SP 227 vented spike, luer lock

Article number

31 9095-100

Filling volume: 0.4 ml

Packing unit: 50 pieces/carton
200 pieces/carton- 11,200 pieces/pallet



Minispike green

31 9099-100

Filling volume: 0.2 ml

Packing unit: 50 pieces/carton
200 pieces/carton- 11,200 pieces/pallet