








INVIVO MR-Monitoring Expression MR400 Accessory Guide

ECG Accessories

Part Number	Description	Image
989803193751	CV ECG Leads, IEC	
989803193761	Standard ECG Leads, IEC	
989803193771	Neonatal ECG Leads, IEC	
989803152291 (9009)	Prep Gel (3 piece/Box) Must be used when applying MRI electrodes	
989803179031 (9303A)	Quadrode, Standard (25 piece/Box)	
989803179041 (9371A)	Quadrode, CV (25 piece/Box)	
989803179051 (9372A)	Quadrode, Neonatal (25 piece/Box)	

SpO2 Accessories










Part Number	Description	Image
989803161991	Quick Connect SpO2 Sensor, MRI	
989803166551	Adult Quick Connect Disposable SpO2 Applicator (20 piece/Box)	
989803166561	Pediatric Quick Connect Disposable SpO2 Applicator (20 piece/Box)	
989803166571	Infant Quick Connect Disposable SpO2 Applicator (20 piece/Box)	
989803166581	Neonatal Quick Connect Disposable SpO2 Applicator (20 piece/Box)	
989803166531	Adult Quick Connect Reusable SpO2 Clip	
989803166541	Pediatric Quick Connect Reusable SpO2 Clip	

NIBP Accessories

Part Number	Description	Image
989803183221	NiBP Hose, Standard (Single Lumen) For use with Adult, Pediatric and Infant Expression MR NiBP Cuffs (excl. Infant #5)	
989803183231	NiBP Hose, Neonatal (Single Lumen) For use with Neonatal and Infant #5 Expression MR NiBP Cuffs	
989803182681	NiBP Cuff, Thigh (44 – 56 cm)	
989803182671	NiBP Cuff, Large Adult-Long (35 – 45 cm)	
989803182661	NiBP Cuff, Large Adult (35 - 45 cm)	
989803182651	NiBP Cuff, Adult-Long (27.5 – 36 cm)	
989803182641	NiBP Cuff, Adult (27.5 – 36 cm)	
989803182631	NiBP Cuff, Small Adult (20.5 – 28.5 cm)	


989803182621	MR200/400 NiBP Cuff, Pediatric (14.0 – 21.5 cm)	
989803182611	MR200/400 NiBP Cuff, Infant (10.0 – 15.0 cm)	
989803183201	Disposable NiBP Cuff, Neonatal #4 (Box of 10) 7.1 – 13.1 cm Delivery time: 2 weeks	
989803183191	Disposable NiBP Cuff, Neonatal #3 (Box of 10), 5.8 – 10.9 cm Delivery time: 2 weeks	
989803183181	Disposable NiBP Cuff, Neonatal #2 (Box of 10), 4.3 - 8 cm Delivery time: 2 weeks	
989803183171	Disposable NiBP Cuff, Neonatal #1 (Box of 10), 3.1 – 5.7 cm Delivery time: 2 weeks	

CO2 and Respiration Accessories (for monitors configured without Anesthetic Agents)


Part Number	Description	Image
989803183241	Disposable LoFlo Sample Line with Adult Cannula (Box of 20)	
989803183251	Disposable LoFlo Sample Line with Pediatric Cannula (Box of 20)	
989803183261	Disposable LoFlo Sample Line with Neonatal Cannula (Box of 20)	
989803183291	Disposable LoFlo Sample Line with Airway Adapter (Box of 20) Delivery time: 2 weeks	
989803185381	Disposable LoFlo Sample Line with Airway Adapter (Box of 100) Delivery time: 2 weeks	
989803183271	Disposable LoFlo Sample Line with Adult Cannula, Divided (Box of 20) Provides inlet for O2 delivery	
989803185361	Disposable LoFlo Sample Line with Adult Cannula, Divided (Box of 100) Provides inlet for O2 delivery	
989803183281	Disposable LoFlo Sample Line with Pediatric Cannula, Divided (Box of 20) Provides inlet for O2 delivery to patient	
989803185371	Disposable LoFlo Sample Line with Pediatric Cannula, Divided (Box of 100) Provides inlet for O2 delivery	

**Anesthetic Agents/O2/CO2 Accessories
(for monitors configured with Anesthetic Agents)**




Part Number	Description	Image
989803152661 (94018)	Anesthetic Agents/O2/CO2 Sample Line	
989803152671 (94012)	Anesthetic Agents/O2/CO2 Water Trap	
989803152561 (9012)	Adult Disposable Cannula For use with 989803152661 Anesthetic Agents/O2/CO2 Sample Line.	
989803152571 (9013)	Pediatric Disposable Cannula For use with 989803152661 Anesthetic Agents/O2/CO2 Sample Line.	
989803152591 (9015)	Infant Disposable Cannula For use with 989803152661 Anesthetic Agents/O2/CO2 Sample Line.	
989803152581 (9014)	Small Infant Disposable Cannula For use with 989803152661 Anesthetic Agents/O2/CO2 Sample Line.	
989803152601 (9016)	Divided Adult Disposable Cannula Provides inlet for O2 delivery. For use with Anesthetic Agents/O2/CO2 Sample Line.	
989803152631 (9016C)	Divided Pediatric Disposable Cannula Provides inlet for O2 delivery. For use with Anesthetic Agents/O2/CO2 Sample Line.	
989803152621 (9016B)	Divided Infant Disposable Cannula Provides inlet for O2 delivery. For use with Anesthetic Agents/O2/CO2 Sample Line.	

<p>989803152611</p> <p>(9016A)</p>	<p>Divided Small Infant Disposable Cannula</p> <p>Provides inlet for O2 delivery. For use with Anesthetic Agents/O2/CO2 Sample Line.</p>	
------------------------------------	--	---


Respiration



Part Number	Description	Image
989803152791	<p>Chest Bellows</p> <p>Connects to port on the Wireless SpO2 Module for measurement of pneumatic respiration.</p>	

Temperature



Part Number	Description	Image
989803194511	<p>FlexTEMP II Temperature Sensor</p> <p>Direct mode axillary, esophageal, or rectal measurement. For esophageal/rectal Measurement, or whenever risk of cross-contamination is a concern, use with 989803172181 FlexTEMP Sensor Jacket.</p>	
989803168891	<p>Temp Sensor Lubricant, 4.25 Oz, (Box of 12)</p> <p>For use during internal applications such as esophageal or rectal measurements.</p>	
989803178181	<p>FlexTEMP Sensor Jacket (Box of 10)</p> <p>For use with 989803162531 FlexTEMP Temperature Sensor or 989803194511 FlexTEMP II Temperature Sensor.</p>	

Invasive Blood Pressure

Part Number	Description	Image
989803194601	<p>Expression IBP Transducer Cable, 5ft</p>	

989803194631	Expression MR IBP Disposable Pressure Transducer Kit, Adult/Pediatric (Box of 20) <i>For use with 989803194601 Expression MR Transducer Cable.</i>	
989803194641	Expression MR IBP Disposable Pressure Transducer Kit, Infant/Neonatal (Box of 20) <i>For use with 989803194601 Expression MR Transducer Cable.</i>	

Power Accessories

Part Number	Description	Image
989803152881	Wireless Module Batteries (3.7V batteries)	
989803191031	Wireless Module Batteries Charger (for charging the 989803152881 batteries)	
989803169491	Monitor Batteries (14.8V batteries)	